

# ELUO™ NEBULIZER CLEANING TOOL

As easy as 1, 2, 3

**1. Fill with solvent**



**2. Insert nebulizer**



**3. Clean nebulizer**



Particle build-up in a nebulizer capillary and tip causes sample flow to be constricted, thereby reducing nebulizer efficiency and performance. Now, blocked nebulizers can be safely and easily restored to optimum performance with the Eluo. The Eluo is designed to efficiently deliver a cleaning solution through the nebulizer capillary to dislodge particle build-up and thoroughly clean the nebulizer. One simple action does it all. No more messy procedures or shattered nebulizers in ultrasonic baths.

Use the Eluo regularly to maintain nebulizer performance and prolong nebulizer life. Every lab should have an Eluo.

Order P/N 70-ELUO for glass concentric nebulizers or P/N 70-ELUO-OPD for OpalMist, PolyCon or DuraMist nebulizers.



P/N 70-ELUO-OPD

## INTERNATIONAL

Glass Expansion • 6 Central Boulevard • Port Melbourne • Vic 3207, Australia  
Telephone: +61 3 9320 1111 • Toll Free (Aust): 1800 777 638  
Facsimile: +61 3 9320 1112 • Email: [enquiries@geicp.com](mailto:enquiries@geicp.com)

## AMERICAS

Glass Expansion • 4 Barlows Landing Road • Unit 2A • Pocasset, MA 02559, USA  
Telephone: 508 563 1800 • Toll Free (US): 800 208 0097  
Facsimile: 508 563 1802 • Email: [geusa@geicp.com](mailto:geusa@geicp.com)

[www.geicp.com](http://www.geicp.com)



**GLASS EXPANSION**  
Quality By Design