

CONEGUARD™ THREAD PROTECTOR


Extend the life of your ICP-MS cones



When cleaning cones which have a screw thread, it is important that the thread does not come in contact with any corrosive solution. If the thread gets corroded, the cone may not seal correctly or it may bond to the base and be difficult to remove. And with Pt cones, the thread is likely to wear out before the Pt insert.

It is also important to keep the thread in good condition to prevent cross-threading and potential damage to the cone housing.

The Glass Expansion ConeGuard™ Thread Protector seals the thread and protects it from corrosion during the cleaning process. The ConeGuard is available for a range of ICP-MS models. Please visit our website or contact us for details.



INTERNATIONAL

Glass Expansion • 6 Central Boulevard • Port Melbourne • Vic 3207, Australia
Telephone: +61 3 9320 1111 • Toll Free (Aust): 1800 777 638
Facsimile: +61 3 9320 1112 • Email: enquiries@geicp.com

AMERICAS

Glass Expansion • 4 Barlows Landing Road • Unit 2A • Pocasset, MA 02559, USA
Telephone: 508 563 1800 • Toll Free (US): 800 208 0097
Facsimile: 508 563 1802 • Email: geusa@geicp.com

www.geicp.com



GLASS EXPANSION
Quality By Design