

Semi Demountable

Glass Expansion Semi Demountable Torches



- Highly accurate construction
- Lowest running cost than 1-piece torches

Usually comprises a quartz torch body, a torch adaptor, an injector and GazFit connectors. All these parts are replaceable, making the semi demountable torch more cost effective than the fixed quartz torch. Essentially the design of the torch for standard analyses and HF analyses is the same, with the exception of the injector material. The torch for aqueous and organics analyses comes with a quartz injector, while the torch version for HF analyses comes with an alumina injector. Hence simply interchange between quartz and alumina injectors for appropriate analyses. Several internal diameters for both quartz and alumina injectors are also available.